



## Thermal analysis of the sample# P9588-AzoMA

Thermal analysis of the samples was carried out on a TA Q100 differential scanning calorimeter at a heating rate of 10°C/min. The midpoint of the slope change of the heat flow plot of the second heating scan was considered as the glass transition temperature ( $T_g$ ).


## Melting and crystallization curve for the sample

The melting temperature ( $T_m$ ) was taken as the maximum of the endothermic peak whereas the crystallization temperature ( $T_c$ ) was considered as the minimum of the exothermic peak.

## Thermal analysis results at a glance

Sample	$T_m$ (°C)	$T_c$ (°C)	$T_g$ (°C)
AzoMA	71/118	66/113	-

## Crystallization curve for AzoMA



## Melting curves for AzoMA

