


Sample Name:


Triglycerol

Sample #: **P9625A-Triglycerol**



Mn x 10 ³ (total)	PDI
0.24	1.10

HNMR of the Triglycerol:



Characterization.

By Size exclusion chromatography (SEC): Varian liquid chromatograph equipped with UV and refractive detector. SEC columns from Supelco were used with THF containing 2 vol% (Et)₃N as the eluent. The molecular weights were determined using light scattering detector and viscosity detector. The molecular weights and the polydispersity indice were calculated.