


Sample Name: Poly(ethylene oxide)-*b*-poly(D,L-lactide)

Sample #: P40627-EOLA (DL form)




#### Composition:

$M_n \times 10^3$ (g/mol) [PEO- <i>b</i> -PLA]	$M_w/M_n$
12.0- <i>b</i> -8.5	1.15

Glass transition temperature, $T_g$ (PEO block):	-42 °C
Cold crystallization temperature, $T_{c.c.}$ :	-17 °C
Melting point, $T_m$ (PEO block)*:	57 °C
* $T_g$ of PLA block overlaps with $T_m$ of PEO block.	

#### Synthesis procedure:

Scheme of poly(ethylene oxide-*b*-lactide) synthesis is shown below:



#### Characterization:

To determine the molecular weight of the first block, an aliquot of anionic poly(ethylene oxide) block was terminated before addition the lactide monomer, and analyzed by size exclusion chromatography (SEC) using DMF as an eluent. The final block copolymer composition was calculated from  $^1\text{H-NMR}$  spectroscopy by comparing the peak area of the methoxy-protons of poly(ethylene oxide) at *ca.* 3.6 ppm and the poly(lactide) protons at *ca.* 5.1 and 1.55 ppm. The polydispersity index ( $M_w/M_n$ ) of the diblock copolymer was obtained by SEC.


#### Thermal analysis:

Thermal analysis was performed on TA Instruments Q100 differential scanning calorimeter (DSC) under a nitrogen atmosphere. The glass transition temperature ( $T_g$ ) and melting point ( $T_m$ ) of the polymer were measured at a scan rate of 10°C/min shortly after creating thermal history of the sample.

#### Solubility:

The polymer is soluble in chloroform, THF, DMF, toluene; and it precipitates from ethanol, ether, and hexanes.


#### $^1\text{H-NMR}$ (500 Mhz, $\text{CDCl}_3$ ) spectrum of PEO-PLA:



#### SEC elugram of the PEO (first block):

P40617-EGOCH3

Conc	33.3957
dndc	0.0440
Solvent	DMF w 0.023M LiBr
Flow Rate	0.7000
Method	PS80k-May2017-0000.vcm




Sample	$M_n$	$M_w$	$M_p$	$M_w/M_n$	IV
P40617_01(553).vdt	11,907	12,349	11,591	1.037	0.1066

#### SEC elugram of the PEO-PLA diblock copolymer:

P40627-EGLA (DL)

Conc	77.7035
dndc	0.0310
Solvent	DMF w 0.023M LiBr
Flow Rate	0.7000
Method	PS80k-May2017-0000.vcm




Sample	$M_n$	$M_w$	$M_p$	$M_w/M_n$	IV
P40627_01.vdt	20,629	23,794	20,053	1.153	0.1065

**DSC thermograms of the polymer  
(2<sup>nd</sup> cooling scan [top] and 2<sup>nd</sup> heating scan [bottom], 10°C/min):**

Sample: P40627-EOLA (DL-form)  
Size: 10.4000 mg

File: P40627\_EO-LA.001



**DSC thermograms of the polymer (3<sup>rd</sup> heating scan, 10°C/min):**

Sample: P40627-EOLA (DL-form)  
Size: 10.4000 mg

File: P40627\_EO-LA.001

