


**Sample Name:** Polystyrene –Broad distribution

**Sample #:** P41447-S

**Structure:**




**Composition:**

Mn x 10 <sup>3</sup>	PDI
294.0	1.5
Tg	100 °C

**Synthesis Procedure:**

Polystyrene is obtained by living anionic polymerization of styrene as illustrated below:



**Characterization:**


The molecular weight and polydispersity index (PDI) are obtained by size exclusion chromatography (SEC) in THF. SEC analysis was performed on a Varian liquid chromatography equipped with refractive and UV light scattering detectors. Three SEC columns from Supelco (G6000-4000-2000 HXL) were used with triple detectors from Viscotek Co.

**Solubility:**

Polystyrene is soluble in DMF, THF, toluene and CHCl<sub>3</sub>. It precipitates from methanol, ethanol, water and hexanes.

**SEC elugram of the polymer in THF:**

Agilent GPC/SEC Software



**DSC thermogram of the polymer:**

