

## Product Profile


### Identification

**Product Name:** Poly(Acrylamide)

**Product Lot Number:** P41550D-AMD

**CAS #:** 9003-05-8

**Chemical Architecture:**




**Composition:**

<b>Mn (g/mole)</b>	<b>6,350,000</b>
<b>Mw (g/mole)</b>	<b>9,463,000</b>
<b>Mw/Mn</b>	<b>1.49</b>
<b>Tg (°C)</b>	<b>184</b>
<b>dn/dc (mL/g) in THF at 30 °C</b>	<b>0.180</b>

### Method of Synthesis

Poly(acrylamide) is synthesized by RAFT polymerization of acrylamide using 4,4'-azo(4-cyanopentanoic acid) as initiator and trithiocarbonate as chain transfer agent in water. The reaction scheme is shown below:



**Solubility in different solvents:**

Water	√
THF	X
Alcohol	X

## Validation of Architecture

### A. Gel Permeation Chromatography (GPC), SEC Profile:

The polymer was characterized by size exclusion chromatography (SEC) using State-of-the-art Agilent Technologies 1260 Infinity II GPC system Equipped with triple detector:

**Solvent (mobile phase)** 2% acetic acid in Millipore water

**Filtration:** 0.45 µNylon Syringe Filter

**Columns:** Agilent three columns


**Flow Rate:** 1 ml/min

**Injection Volume:** 100 µL

**Column Temperature:** 30 °C


**Calibration of Instrument using PEO polymer.**

**Note:** Polyacrylamide bearing Mw > 1M are difficult to filter therefore this equipment is highly sensitive where less than 1mg/ml polymer solution can be detected by triple detector.



MW Extrapolated Fit Calculated Statistics (g/mol)			
MW	9463248	Mn	6350725
Mp	4771044	Mz	21076480
Mz+1	41510562	Mv	16040008
PD	1.49		
Structural Statistics (Rg)			
Rgw	113.82	Rgn	109.61
Rgp	106.69	Rgz	122.75
Rgz+1	130.85		

### B. DSC thermogram for the polymer:



124 avenue Avro, Dorval (Montreal)  
 Quebec H9P 2X8 Canada  
 Phone : +1-514-421-5517 or +1-514-421-5506  
[support@polymersource.com](mailto:support@polymersource.com)  
[www.polymersource.ca](http://www.polymersource.ca)