


Sample Name: Polystyrene

Sample #: P4555-S

Structure:




Composition:

$M_n \times 10^3$	PDI
217.0	1.09

Synthesis Procedure:

Polystyrene is obtained by living anionic polymerization of styrene as illustrated below:



Characterization:


The molecular weight and polydispersity index (PDI) are obtained by size exclusion chromatography (SEC) in THF. SEC analysis was performed on a Varian liquid chromatograph equipped with refractive and UV light scattering detectors. Three SEC columns from Supelco (G6000-4000-2000 HXL) were used with triple detectors from Viscotek Co.

Solubility:

Polystyrene is soluble in DMF, THF, toluene and CHCl_3 . It precipitates from methanol, ethanol, water, and hexanes.

SEC profile of the Sample:

P4555-S



Size Exclusion Chromatography of polystyrene

$M_n=217000$, $M_w=237000$, $M_w/M_n=1.09$

Thermal analysis of the sample was carried out using a differential scanning calorimeter (TA Q100) at a heating rate of 10°C/min . The inflection glass transition temperature (T_g) has been considered.

