


Sample Name: Poly(Vinyl alcohol)

SEC of Homopolymer:

Sample #: P11285-VA

P11285-VA

Structure:




Composition:

$M_n \times 10^3$	PDI
110.0	1.8


Synthesis Procedure:

Poly vinyl alcohol is obtained from the hydrolysis of poly vinyl acetate under alkaline conditions. Poly vinyl acetate is obtained by free radical polymerization using iodo ethyl acetate as chain transfer reagent. Polymerization was carried out in bulk. The polymerization scheme can be illustrated as follows:



Characterization:

The molecular weight and polydispersity index (PDI) are obtained by size exclusion chromatography (SEC) in THF. SEC analysis was performed on a Varian liquid chromatograph equipped with refractive and UV light scattering detectors. Three SEC columns from Supelco (G6000-4000-2000 HXL) were used with triple detectors from Viscotek Co. ^1H NMR analysis was carried out on Varian instrument at 500MHz.



Size exclusion chromatograph of Poly vinyl acetate:

$M_n=215,000$ $M_w=388,000$ $PI=1.8$

After Hydrolysis of acetate group: M_n 110,000 M_w/M_n : 1.8