


Sample Name: Poly(methyl methacrylate)

Sample #: P20205-MMA

Structure:



Composition:

Mn x 10 ³	PDI
145.0	1.55
Syndio : Hetero : Iso	65:32:3
T _g	119°C

Synthesis Procedure:

Tacticity of the poly(methyl methacrylate) is tailored by anionic polymerization of MMA monomer in different polarity solvents mixture and using different ligands.

Characterization:


Tacticity of the polymer was determined by ¹H NMR. The molecular weight and polydispersity index (PDI) were obtained by size exclusion chromatography (SEC) in THF. SEC analysis was performed on a Varian liquid chromatograph equipped with refractive and UV light scattering detectors. Three SEC columns from Supelco (G6000-4000-2000 HXL) were used with triple detectors from Viscotek Co.

Thermal analysis of the samples was carried out using a differential scanning calorimeter (TA Q100) at a heating rate of 10°C/min. The inflection glass transition temperature (T_g) of the sample has been considered.

Solubility:


The polymer is soluble in chloroform.

¹H NMR spectrum of PMMA:




SEC elugram of PMMA homopolymer: Sample ID:P20205-MMA

Concentration (mg/mL)	8.0490
Sample dv/dc (mL/g)	0.0840
Method File	PS80K-March6-2015-0000.vsm
Column Set	3x PL 1113-6300
Solvent	THF




Sample	MW Number Average (Da)	MW Weight Average (Da)	MW at Peak (Da)	Polydispersity	Intrinsic Viscosity (dL/g)
P20205-MMA_01(517).vdt	145,284	223,586	216,288	1.539	0.3734

DSC thermogram :




Thermograms of PMMA:


(a) syndiotactic >79%




(d) atactic



(b) syndiotactic >85%



(c) isotactic >97%



Summary of DSC results for PMMA of different tacticity:

PMMA microstructure	Tacticity Syndio : Iso : Hetero	T_g (°C)
Syndiotactic >79%	79 : 19 : 2	120
Syndiotactic >85%	86 : 0 : 14	123
Isotactic >97%	0 : 97 : 3	44
Atactic	56 : 6 : 38	97