




¹H-NMR Spectrum of the random copolymer before hydrolysis:



SEC of the random copolymer before hydrolysis:



DSC thermogram for the sample:




Sample Name:

Random Copolymer Poly(acrylamide-co-N,N-dimethyl acrylamide)

Sample #: P4398-AMDNNDMA

Structure:



Composition:

PAMD (mol%) : 56

$M_n \times 10^3$	PDI
1000.0	5.0
T_g for the random copolymer	189°C

Synthesis Procedure:

The polymer is prepared by RAFT process.

Characterization:

The molecular weight and polydispersity index (PDI) were calculated by SEC using water as eluent. The copolymer composition was calculated from ¹H-NMR spectroscopy.

Thermal analysis:

Thermal analysis of the samples was carried out on a TA Q100 differential scanning calorimeter at a heating rate of 20°C/min. The midpoint of the slope change of the heat flow plot of the second heating scan was considered as the glass transition temperature (T_g).

Solubility:

The polymer is soluble in acetone, insoluble in ether, hexane.