


Sample Name: Poly(ethylene oxide)-*b*-poly(D,L-lactide)

Sample #: P40613-EOLA (DL form)




Composition:

$M_n \times 10^3$ (g/mol) [PEO- <i>b</i> -PLA]	M_w/M_n
11.0- <i>b</i> -15.0	1.50

Glass transition temperature, T_g (PEO block):	-16 °C
Cold crystallization temperature, $T_{c.cr.}$:	32 °C
Melting point, T_m (PEO block)*:	56 °C
* T_g of PLA block overlaps with T_m of PEO block.	

Synthesis procedure:

Scheme of poly(ethylene oxide-*b*-lactide) synthesis is shown below:



Characterization:

To determine the molecular weight of the first block, an aliquot of anionic poly(ethylene oxide) block was terminated before addition the lactide monomer, and analyzed by size exclusion chromatography (SEC) using DMF as an eluent. The final block copolymer composition was calculated from $^1\text{H-NMR}$ spectroscopy by comparing the peak area of the methoxy-protons of poly(ethylene oxide) at *ca.* 3.6 ppm and the poly(lactide) protons at *ca.* 5.1 and 1.55 ppm. The polydispersity index (M_w/M_n) of the diblock copolymer was obtained by SEC.


Thermal analysis:

Thermal analysis was performed on TA Instruments Q100 differential scanning calorimeter (DSC) under a nitrogen atmosphere. The glass transition temperature (T_g) and melting point (T_m) of the polymer were measured at a scan rate of 10°C/min shortly after creating thermal history of the sample.

Solubility:


The polymer is soluble in chloroform, THF, DMF, toluene; and it precipitates from ethanol, ether, and hexanes.

$^1\text{H-NMR}$ (500 Mhz, CDCl_3) spectrum of PEO-PLA:



SEC elugram of the PEO (first block):

P5675-EGOCH3




Size exclusion chromatograph of Poly(ethylene glycol) methyl ether:
 $M_n=11000$, $M_w=12000$, $PI=1.09$

SEC elugram of the PEO-PLA diblock copolymer:

P40613-EOLA

Conc	34.4029
dn/dc	0.0330
Solvent	DMF w 0.023M LiBr
Flow Rate	0.7000
Method	PS80k-Mey2017-0000.vcm




Sample	M_n	M_w	M_p	M_w/M_n	IV
P40613-EOLA_01(533).vdt	25,222	37,859	25,305	1.501	0.1583

DSC thermograms of the polymer (2nd heating scan, 10°C/min):

[Note: no crystallization peak is observed for cooling scans.]

Sample: P40613_EOLA

File: D:\00\P40613_EO-LA.001
Operator: Iryna
Run Date: 29-May-2017 09:04



Sample: P40613_EOLA

File: D:\00\P40613_EO-LA.001
Operator: Iryna
Run Date: 29-May-2017 09:04

